


# HANAR STENT®

## TwynLay™

COLON/RECTUM  
DUODENUM/PYLORUS





### “CONFORMS ANYWHERE”

**For treatment of malignant obstruction  
in Colon/Rectum and Duodenum/Pylorus**




# HANAR STENT®

## TwynLay™

COLON/RECTUM  
DUODENUM/PYLORUS



- ☑ Silicone membrane in middle of inner and outer mesh
- ☑ Silicone membrane blocks tissue in-growth
- ☑ Outer mesh minimize the risk of migration



## Ordering Information

※ Sizing and availability varies by country

### • TLC for Colon/Rectum | Endoscopic Application

| Model          | Stent (mm) |                |               | Delivery Device |                  |
|----------------|------------|----------------|---------------|-----------------|------------------|
|                | Diameter   | Usable Length* | Total Length* | Length (mm)     | Diameter (mm/fr) |
| TLC-20-060-230 | 26-20-26   | 30             | 60            | 2300            | 3.4/10.2         |
| TLC-20-150-230 |            | 120            | 150           | 2300            | 3.4/10.2         |
| TLC-22-060-230 | 28-22-28   | 30             | 60            | 2300            | 3.4/10.2         |
| TLC-22-150-230 |            | 120            | 150           | 2300            | 3.4/10.2         |

### • TLC for Colon/Rectum | Fluoroscopic Application

| Model          | Stent (mm) |                |               | Delivery Device |                  |
|----------------|------------|----------------|---------------|-----------------|------------------|
|                | Diameter   | Usable Length* | Total Length* | Length (mm)     | Diameter (mm/fr) |
| TLC-20-060-140 | 26-20-26   | 30             | 60            | 1400            | 3.4/10.2         |
| TLC-20-150-140 |            | 120            | 150           | 1400            | 3.4/10.2         |
| TLC-22-060-140 | 28-22-28   | 30             | 60            | 1400            | 3.4/10.2         |
| TLC-22-150-140 |            | 120            | 150           | 1400            | 3.4/10.2         |

### • TLD for Duodenum/Pylorus | Endoscopic Application

| Model          | Stent (mm) |                |               | Delivery Device |                  |
|----------------|------------|----------------|---------------|-----------------|------------------|
|                | Diameter   | Usable Length* | Total Length* | Length (mm)     | Diameter (mm/fr) |
| TLD-20-060-230 | 26-20-26   | 30             | 60            | 2300            | 3.4/10.2         |
| TLD-20-150-230 |            | 120            | 150           | 2300            | 3.4/10.2         |
| TLD-22-060-230 | 28-22-28   | 30             | 60            | 2300            | 3.4/10.2         |
| TLD-22-150-230 |            | 120            | 150           | 2300            | 3.4/10.2         |

### • TLD for Duodenum/Pylorus | Fluoroscopic Application

| Model          | Stent (mm) |                |               | Delivery Device |                  |
|----------------|------------|----------------|---------------|-----------------|------------------|
|                | Diameter   | Usable Length* | Total Length* | Length (mm)     | Diameter (mm/fr) |
| TLD-20-060-120 | 26-20-26   | 30             | 60            | 1200            | 3.4/10.2         |
| TLD-20-150-120 |            | 120            | 150           | 1200            | 3.4/10.2         |
| TLD-22-060-120 | 28-22-28   | 30             | 60            | 1200            | 3.4/10.2         |
| TLD-22-150-120 |            | 120            | 150           | 1200            | 3.4/10.2         |

\* Increment of 10mm