





Reloadable Hemoclip System

- Easy & rapid clip exchange
- Various specifications
- 1:1 rotation
- Re-openable clip




Indication

This clip is compatible with endoscope, which is indicated for clip placement within digestive tract for the purpose of mechanical pressure treatment of bleeding of small arteries and pulsation. The clip is intended for single use.



Easy & rapid clip exchange

A simple push and pull, saving time and effort, significantly improving surgical efficiency and cost-effectiveness




Various specifications

Available in multiple opening widths, meeting different clinical needs



1:1 rotation

Bidirectional synchronous rotation for precise control and easy operation



Re-openable clip

Easy repositioning ensuring accurate placement

Ordering Information

Product	Product Code	Cartridge Color	Opening Span	Claw Angle	Clip O.D.	Working Length	Rotatable	Repositionable	Package Unit
EndoAll Hemoclip Applicator	AG-55501-2300-135-11	-	11mm	135°	2.6mm	2,300mm	Y	Y	10
EndoAll Hemoclip Applicator	AG-55501-2300-135-16	-	16mm	135°	2.6mm	2,300mm	Y	Y	10
EndoAll Hemoclip Reload	AG-55503-0000-135-11		11mm	135°	2.6mm	-	-	Y	25
EndoAll Hemoclip Reload	AG-55503-0000-135-16		16mm	135°	2.6mm	-	-	Y	25