


EndoAll Guidewire

- Nitinol core construction
- Hydrophilic tip
- Dual-color spiral pattern






Nitinol core construction

Optimized stiffness for kink resistance, easier access, and precise control

Hydrophilic tip



Ultra-smooth, highly flexible tip for safe and easy passage through tight stenosis



Dual-color spiral pattern

Smooth tracking and easy device exchange with reliable guidance

Ordering Information

Product Image	Product Code	Tip Type	Outer Diameter	Working Length	Hydrophilic tip coating length
	AG-5043-2545	Straight	0.025inch (0.63mm)	4,500mm	50mm
	AG-5043-3545	Straight	0.035inch (0.88mm)	4,500mm	50mm
	AG-5044-2545	Angled	0.025inch (0.63mm)	4,500mm	50mm
	AG-5044-3545	Angled	0.035inch (0.88mm)	4,500mm	50mm